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
Data Acquisition Workstation	
Operating System & Version	Microsoft Windows™ 7 Professional with Service Pack 1 (64-bit) English
System Software	Thermo Scientific Xcalibur™ v4.1 Thermo Scientific FreeStyle™ 1.3 Thermo Scientific mzVault™ 1.0 SP1 Thermo Scientific TraceFinder™ Core 4.1 SP2 Thermo Scientific TraceFinder™ General Quan 4.1 SP2 Shimadzu Instrument driver 6.5 Thermo Scientific Foundation™ v3.1 Thermo Scientific Exactive Series™ 2.9 Microsoft.NET Framework 4.7 Microsoft Office Pro Plus 2016

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Data Acquisition Workstation	
Operating System & Version	Microsoft Windows™ 10 Professional (64-bit) English
System Software	Xcalibur v4.1 Foundation v3.1 Free Style v1.3 TraceFinder Core 4.1 SP2 TraceFinder Generral Quan 4.1 SP2 Thermo Scientific Compound Discoverer™ v2.1 Microsoft Office Standard 2016

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Windows Folder	Group Name			
	Xcal_ADMIN	Xcal_SCI	Xcal_QA	Xcal_PM
C:\Xcalibur	Full Control	Read/Write/ExecuteList	Read only	Not configured
C:\Xcalibur\data		Folder Contents		
C:\Xcalibur\Methods		Read/Write/Execute /Modify		
C:\Xcalibur\PM_Service		Read/Write /Execute List Folder Contents		
C:\Program Files\Thermo\Foundation				
C:\Program Files(x86)\Thermo\Foundation\AuditSystemLogs				
C:\Program Data\Thermo Scientific\INI				

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5 → Test Instructions


5.1 → ISO 2, SAqX 1, SAqX 2

Initiate and Complete a Sample Acquisition Sequence

Confirm that the system can initiate a sequence sample acquisition to allow only authorized users to securely acquire data to a local storage location.

1. → Log onto the workstation with an account in the **Xcal SCI** role.

2. → Open the **Tune**, then set the following figure.

		actual
Sheath gas flow rate	50	0
Aux gas flow rate	10	0
Sweep gas flow rate	0	0
Spray voltage (kV)	3.50	0.00
Spray current (µA)		0.00
Capillary temp. (°C)	325	325
S-lens RF level	50.0	
Aux gas heater temp (°C)	350	36

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